

ENTFirst

Enables clear observation of nystagmus

Frenzel Goggles S5

Japanese technology and spirit cultivated for
over 60 years, distilled into
one exceptionally reliable Frenzel Goggles

Made in Japan



Enables clear observation of nystagmus

Frenzel Goggles S5



Bulb-colored LEDs on right and left



Clear observation for eye movement



Comfort

- Well-fitting design for patient's face and operator's hands
- Bulb-colored LEDs reveal the eye movement in familiar color



Easy to Use

- Lightweight and easy to carry (only 170 g without fixing belt)
- Easy to replace face pad
- Single AAA battery operation



Reliable

- High quality lenses allow clear observation of nystagmus
- Durable and lightweight ABS resin frame allows easy handling

Specifications

Lenses	High quality convex plastic
Illumination	2 bulb-colored LEDs (right and left side)
Frame	ABS resin
Fixing belt	Artificial leather velcro
Face pad	Replaceable
Power source	One AAA battery
Weight	195 g (with fixing belt) 170 g (without fixing belt)
Rated voltage	1.5V DC
Current consumption	0.17 A

ENTFirst

DAIICHI MEDICAL CO., LTD. www.first-med.co.jp/en/

2-27-16, Hongo, Bunkyo-ku, Tokyo 113-0033, Japan Tel: 81-3-3814-0111 Fax: 81-3-3814-0135 Email: first-ex@first-med.co.jp

022012

Color on printed materials differ from color of actual product. Design and specifications are subject to change without prior notice.